



## MINIMAX-VAC

### Features

- » 1 High-performance suction hose 220 l/min
- » 1 rotary knob for continuously variable regulation
- » Patent-pending BPR Swiss technology
- » With adapter for saliva ejector tips and oral-dry-cup
- » 1 Suction fluid container: 1 litre (automatic overflow protection)
- » 1 DENTA-BAG carrying bag

### Technical Data

- » Weight: 3.9 kg (weight including DENTA-BAG 4.5 kg)
- » Dimensions MINIMAX-VAC: 190mm x 140mm x 285mm
- » Dimensions DENTA-BAG carrying bag: 275mm x 300mm x 275mm

**Warranty:** 2 years

**Options:** DENTA-PRO, DENTA-TROLLEY, DENTA-CASE etc.

## HIGHLIGHTS

CONTINUOUS OPERATION

3.9 KG

FULL PERFORMANCE

SUCTION 220 L/MIN

DENTA-BAG

PLUG-AND-PLAY

The MINIMAX-VAC by BPR.<sup>+</sup>  
High-tech in miniature.

The MINIMAX-VAC sets new standards in portable plug-and-play suction units. To start the MINIMAX-VAC it only needs to be plugged into mains power. With a weight of only 3.9 kilograms (!), the MINIMAX-VAC, with over 220 l / min (!), achieves a comparable suction to a permanently installed treatment unit in a dental office. The suction can be seamlessly regulated any time. Due to this pioneering key data MINIMAX-VAC can be used in many areas such as mobile dentistry, in the dental office, or in numerous other areas of medicine.

The MINIMAX-VAC from BPR. High-tech in miniature.

