

Baby Anne™

This low-cost, lifelike infant CPR trainer is the perfect supplement to the Resusci Baby manikin



Marketed: Worldwide

Type: Full-body infant

Airway: Re-breathing, disposable

Carotid Pulse: No

Squeeze Bulb: n/a

Compr. Generated: n/a

Ribs, Xiphoid: Yes

Gastric, distension simulation: No

Navel: No

Tongue/teeth: No

Adam's Apple: No

Chart Recorder: No

Lightbox: No

Metronome: No

Length: 50 cm (20 in)

Weight: 1.5 kg (3.3 lb)

Packed weight: 1.8 kg (4 lb)

Dimensions packed: 50 cm x 15 cm x 20 cm (20 in x 6 in x 8 in)

Warranty: 1 year

Other Specifications: Lightweight with realistic features. Supplement to Resusci Baby. Softpack. Choking kit available. 4-pack version available.