

PHILIPS

RESPIRONICS

DreamStation

Rediscover
dreams



The future of **sleep therapy** is now

Obtaining a good night's sleep is essential for patients to feel refreshed and energetic.

The adoption of collaborative and connected care has made the dream of a good night's sleep more achievable than ever before. However, sleep therapy must also accomplish today's sleep challenges;



Increase patient adherence

Long-term adherence can be predicted from the first days of CPAP use*. So a first positive experience is crucial for therapy adoption.



Decrease operational costs

Remote troubleshooting tools can help reduce home visits, but also the number of normally functioning devices that are returned for service.



Improve resources utilization

Direct connections to Home Sleep Testing devices offer the opportunity to perform polygraphies at the patient's home without requiring the use of the sleep lab's beds.



Deliver patient education

Patient education is key to therapy success. Easy-to-use, comfortable and intuitive therapy solutions can facilitate patient education and help save precious time.

*Aloia and al., Behav Sleep Med 2007, How early in treatment is PAP adherence established? Revisiting night-to-night variability.

Smaller, lighter, quieter¹

Introducing DreamStation, our multifunctional, connected sleep apnea therapy device that adapts to the unique needs of your patients, supporting therapy adherence from the initial adoption phase through to the integration of therapy into your patient's lifestyle.

A positive airway pressure (PAP) sleep therapy device that offers exceptional comfort, intelligent monitoring, adjustment, and flexible connectivity to the care team for proven clinical care and effortless patient management.

DreamStation is less medical looking than its predecessors and patients have demonstrated overwhelming support, with 95% of trial participants indicating they like the appearance².



Almost 20% smaller and 30% lighter than System One³.



Slim profile makes the DreamStation discrete and easy to transport.



Day to day progress tracking, helps long term adherence.

¹ Compared to System One.

² 2015 Philips User Preference trial, data on file.

³ 2015 Philips internal testing, data on file.



“ Very impressed.
It looked like a
dream machine.”

“ Good-looking piece of
equipment! It’s modern
looking, doesn’t look
like a medical machine.
It’s quieter, too.”

Insight driven innovation

It's easy to develop products using the latest technology but true innovation is only achieved when the solution is developed with a real understanding of the users' needs. This is why we've spent the past 5 years researching, developing and challenging what we know about sleep therapy.

- 700 surveys, interviews and user trials undertaken
- In depth research and development completed across 3 continents
- Developed using the help and experience of physicians, providers and patients

What research told us

- The therapy device needs to be easy to use, discreet and shouldn't look like a medical device – patients don't want to be reminded that they are sick every night
- The device needs to be transportable
- CPAP does not feel natural – patients want an easier way to get used to CPAP therapy
- CPAP pressure can be hard to tolerate all night
- Patients want to know how well they did each morning – they also want to troubleshoot their device remotely



That is why we've developed the DreamStation, an effective sleep therapy device, delivering innovative features that encourage patients to take an active role in their sleep therapy. Giving you the peace of mind that your patients are adhering to their therapy and achieving an exceptional night of sleep.



Gentle comfort

Philips' unique pressure acclimation function, 'EZ Start', offers automatic, personalized adjustments to CPAP pressure which can facilitate CPAP therapy adoption. Equipped with the new 'SmartRamp' feature, DreamStation delivers a reduced pressure while the patient is falling asleep, but still automatically detects and resolves respiratory events.



Easy to use

Intuitive navigation with a display that can be operated from bed, offering patients the ability to track their progress every morning. DreamStation humidifier has been designed to make it easy for patients to use daily and is heated tube capable.



A helping hand

The Performance Check function is designed to make troubleshooting easier than ever before. It helps you to quickly determine which devices require repair so you can avoid wasting time on devices that are running correctly.



Ordering information

DreamStation CPAP	Part number
DreamStation CPAP	INX200S15
DreamStation CPAP with humidifier	INX200H15
DreamStation CPAP with humidifier and heated tube	INX200T15
DreamStation CPAP Pro	
DreamStation CPAP Pro with humidifier	INX400H15
DreamStation CPAP Pro with humidifier and heated tube	INX400T15
DreamStation Auto CPAP	
DreamStation Auto CPAP	INX500S15
DreamStation Auto CPAP with humidifier	INX500H15
DreamStation Auto CPAP with humidifier and heated tube	INX500T15
DreamStation BiPAP Pro	
DreamStation BiPAP Pro	INX600S15
DreamStation BiPAP Pro with humidifier	INX600H15
DreamStation BiPAP Pro with humidifier and heated tube	INX600T15
DreamStation Auto BiPAP	
DreamStation Auto BiPAP	INX700S15
DreamStation Auto BiPAP with humidifier	INX700H15
DreamStation Auto BiPAP with humidifier and heated tube	INX700T15

Humidification	Part number
DreamStation humidifier	EUXH
DreamStation water tank	1122520
RP-DS Humidifier Dry Box Assembly	1120668
RP-DS Humidifier Flip Lid Seal	1120617
RP-DS Humidifier Dry Box Inlet Seal	1120613
Power supply	
DreamStation 80-watt power supply	1118499
Power cord	1005894
Tubes	
15MM Std Tube – DreamStation - RP	PR15
15MM Heated Tube	HT15
Filters	
Pollen filter, reusable (1 per pack)	1122446
Ultra-fine filter, disposable (1 per pack)	1122447
Ultra-fine filter, disposable (2 per pack)	1122518
Ultra-fine filter, disposable (6 per pack)	1122519
Oximetry	
Link module	1120293
Nonin SpO ₂ assembly	1121694
Reuseable finger sensor (clip)	936
Oximetry finger sensor, adult, flex	953
Flexiwrap, sens tape, adult, 25/pk	954A
Reuseable finger sensor (soft), medium	1062050
Reuseable finger sensor (soft), small	1062051
Travel accessories	
Universal PAP/laptop travel briefcase	1114784
DreamStation replacement carrying case	1121162

