




# MAKE IT SAFER

Autoclave  
Steam Sterilizer  
Premium+ Line



- ensure device to work safely
- ensure device to work reliably



- protect users to operate safely
- protect instruments from damage

- ensure complete sterilization
- avoid patients from cross infection


# Double Load



- With the particular 5-layer tray or s/s basket at your option

- Not only the instruments can be laid randomly, but more quantity can be loaded

- To make your work easier more efficient and costs reduced



# Super Fast

- It only takes 39 minutes of the fastest cycle with the full load of 60 handpieces.
- Designed for the high efficiency of instruments reprocessing.

LAFOMED

# Fast Dry Autoclave

**TSD** Technology of Seppy Dry



- Equipped with our state-of -the-art Technoloty of Speedy Dry (TSD)
- The instruments are competely dry.

# NO wet instruments

## Premium+ 18L

Chamber Size  
**Φ 25\*35cm**

G.W./ N.W.  
**63 / 55 kgs**

Product Size  
**60\*50\*43cm**

Power  
**50HZ 1.95kW**



## Premium+ 23L

Chamber Size  
**Φ 25\*45cm**

G.W./ N.W.  
**65 / 57 kgs**

Product Size  
**60\*50\*43 cm**

Power  
**50HZ 1.95kW**

### Program

134 °C Universal

134 °C Quick Universal

121 °C Universal

134 °C Prion

Test Vacuum

Test B-D/Helix