

ID 213

Instrument disinfection



Aldehyde-free concentrate for cleaning and disinfecting both general and surgical instruments and rotating instruments that are sensitive to alkali and alcohol in practices, laboratories and clinics in the dental sector.



wide spectrum of action

fast-acting

Contact time

Recommendation



immersion
bath

5 min

2 min

ultrasonic

ideal for use in
ultrasonic devices

economical
(2.5 l concentrate = 125 l
ready-to-use solution)

high corrosion protection
and certified
material compatibility

recommended by leading
instrument manufacturers

Spectrum of action *



High level of protection against infection:

Bactericidal, tuberculocidal, mycobactericidal in an ultrasonic device, levurocidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV), noro-viruses and polyomaviruses SV40



Key features

14 day
service life and low
application concentration

a single
concentrate for all
disinfection and
cleaning processes

optimum material
compatibility
thanks to
special corrosion
inhibitors



*Tested in accordance with EN 13727, EN 13624, EN 14348, EN 14476, EN 14561, EN 14562, EN 14563 and EN 17111.