

ID 220

Bur disinfection



Aldehyde-free, ready-to-use solution for cleaning and disinfecting rotary instruments in practices, laboratories and clinics in the dental sector.



ready-to-use

very fast-acting

exceptional
cleaning effect

high corrosion
protection

tested material
compatibility

Contact time



immersion
bath

1 min

30 sec

ultrasonic

Recommendation

dissolves **even the most stubborn contamination** on drills and cutters made from alkali-resistant steel and stainless steel, carbide

Spectrum of action *



High level of protection against infection:

Bactericidal, mycobactericidal, tuberculocidal, yeasticidal, fungicidal, virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses, as well as non-enveloped viruses such as adenoviruses, noroviruses, polioviruses, polyomaviruses SV40)



Key features

fully
virucidal

exceptional
cleaning power,
particularly for
the microstructures
in drills and cutters

excellent material
compatibility
thanks to
special corrosion
inhibitors



*Tested in accordance with EN 13727, EN 13624, EN 14348, EN 14476, EN 14561, EN 14562, EN 14563 and EN 17111.