

# FD 333

## FD 333 wipes / FD 333 top wipes

### Quick-acting disinfection



Aldehyde-free, alcohol-based ready-to-use solution for quick-acting disinfection of the surfaces of invasive and non-invasive medical devices and fixtures in practices, laboratories and clinics.



tested using the  
4-field test (EN 16615)

fast-acting with  
a high level of protection  
against infection

rapid and residue-free  
drying

recommended by leading  
equipment manufacturers

available as a ready-to-  
use solution and wipes

#### Contact time



only

# 1 min

fast availability of  
treatment chair

#### Recommendation

for quick-acting  
disinfection of the  
surfaces of medical  
devices and fixtures

#### Spectrum of action \*



**High level of protection against infection:**

Bactericidal, tuberculocidal, fungicidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses)



#### Key features

limited virucidal

# PLUS

quick-acting and  
long-lasting  
disinfection and  
cleaning

good material  
compatibility  
with  
lemon freshness



\*Tested in accordance with EN 13727, EN 13624, EN 14348, EN 14476 and EN 16615.