

FD 366 sensitive / FD 366 sensitive wipes / top wipes

Disinfection of sensitive surfaces



Aldehyde-free, alcohol-based ready-to-use solution for quick-acting disinfection of the sensitive surfaces of invasive and non-invasive medical devices in practices, laboratories and clinics.



ready-to-use and
fast-acting

tested using the
4-field test (EN 16615)

minimum alcohol content

tested material
compatibility

available as a ready-to-
use solution and wipes

Contact time

Recommendation



only

1 minute

fast availability of
treatment chair

ideal for
sensitive surfaces
(plastics, acrylic
glass or
imitation leather)

Spectrum of action *



High level of protection against infection:

Bactericidal, tuberculocidal, yeasticidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), plus noroviruses



Key features

only
17%
alcohol content

maximum
material
compatibility

long-term
maintenance
of the value of
surgery inventory



*Tested in accordance with EN 13727, EN 13624, EN 14348, EN 14476 and EN 16615.