

MD 520

Impression disinfection



Highly effective formaldehyde-free ready-to-use solution for cleaning and disinfecting dental impressions, impression trays, prosthetic work etc., for immersion disinfection in the Dürr Dental Hygobox / disinfection box.



fully virucidal

suitable for universal use

alcohol-free

increased precision of plaster master model

preservation of dimensional stability / plaster compatibility

Contact time

Recommendation

 only **5 min** in accordance with VAH

can be used for impression materials, prosthetic work, orthodontic appliances, drilling templates for implant systems, etc.

Spectrum of action *

 **High level of protection against infection:** Bactericidal, tuberculocidal, yeasticidal, virucidal (enveloped viruses incl. HBV, HCV, HIV, coronaviruses, as well as non-enveloped viruses such as adenoviruses, polyomaviruses SV 40, polioviruses, noroviruses)



Key features



*Tested in accordance with EN 13727, EN 13624, EN 14348, EN 14476, EN 14561, EN 14562, EN 14563 and EN 17111.