

Orotol plus

Disinfection of suction systems



Foam-free liquid concentrate for the simultaneous disinfection, deodorisation, cleaning and care of dental suction units, amalgam separators and spittoons.



maintains the value of your suction system

tested material compatibility

recommended by leading equipment manufacturers

maximum effectiveness when used with cold water

works perfectly in combination with MD 555 cleaner

Contact time



2%

5-60 min

once or twice daily

Recommendation

safe **protection against infection** and **effective cleaning** of organic deposits

Spectrum of action*



High level of protection against infection:

Bactericidal, yeasticidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), as well as adenoviruses and noroviruses



Key features

only **2%** application concentration

disinfects and dissolves biofilm

guarantees technically and hygienically perfect operation of the suction system



*Tested in accordance with EN 13727, EN 14561, EN 13624, EN 14562, EN 14476, EN 17111 and EN 16615.